



## Certificate of Quality

### RosetteArray™ Plate

**LOT #: 13.16.12072023**

Product description:

- Clear 96-well
- 250µm diameter circular micropattern array

The following inspections were carried out and passed:

Visual inspection of each micropattern surface and plate upon assembly.

One plate out of the batch, chosen at random, displayed restriction of cell adherence 3 days after cell seeding for 60 out of 60 of the interior wells.

Certification date:

12-12-2023

Certification Signature:

A handwritten signature in blue ink, appearing to be "D. J. ...", written over a horizontal line.