



Certificate of Quality

LOT #: 20251018.51

Product description: RosetteArray plate

- Clear 96-well
- 250 μ m diameter circular micropattern array

The following inspections were carried out and passed:

Visual inspection of each micropattern surface and plate upon assembly.

One plate out of the batch, chosen at random, displayed restriction of cell adherence 7 days after cell seeding for 60 out of 60 of the interior wells.




Fig. 1. Top) BF images taken at low mag. on Day 4 display tissue restriction. Bottom) BF images taken on Day 7 display typical rosette tissue morphology.

Certification date:

10/31/2025

Certification Signature:

gk